



# Medical mask

Disposable, nonsterile, face mask

# TECHNICAL SHEET

## PRODUCT DESCRIPTION:

Medical mask for human use to cover the face (mouth, nose and chin), made of 3 layers of medical grade protection: a non-woven layer, a middle layer - filter obtained by melting and blowing and an inner layer adapted to skin contact. The lightweight elastic straps for the ears allow the mask to be easily put on or removed.

## COMPOSITION:

- A non-woven layer
- A middle layer - filter obtained by melting and blowing
- An inner layer adapted to skin contact
- Elastic straps for the ears.
- Aluminum clips for the nose, which allows you to adjust the mask.


## PERFORMANCE:

- Bacterial Filtration Efficiency
- Breathable
- Projection resistant.

## FEATURES:

- Bacterial Filtration Efficiency (BFE) (%)  $\geq 95\%$
- Differential pressure (Pa/cm<sup>2</sup>) < 40
- Projection resistance pressure (KPa) Unsked
- Microbial property (cfu/g)  $\leq 30$
- Non-sterile

## OPERATING INSTRUCTIONS:





## DIMENSIONS:

17,5 x 9,5 mm - mask for adults

## PACKAGING:

Box of 50 masks



  Bernard  
comanda@vvup.org  
+40 738 886 499

  Mariana  
comanda@vvup.org  
+40 760 665 659

 Ray  
comanda@vvup.org

# FROM THE MANUFACTURER TO THE CUSTOMER