

# V-VEIL UP<sup>2</sup> SELF SAMPLING VEIL COLLECTOR

## your solution for HPV & STI detection



### SELF-SAMPLING VEIL ADVANTAGES

- ✓ Specially adapted to women who feel discomfort during gynecological examinations, which may also be related to cultural or religious values.
- ✓ The Patented Veil Technology is able to release 90.7% of DNA and 81.3% of Proteins.  
Other methods of self-sampling are extremely FAR of this efficiency.
- ✓ 97.1% acceptance rate on a panel of 253 women.

### IF MORE DNA IS RELEASED, YOU SAVE MORE LIVES!

- ✓ The use of this veil makes it possible to almost double the detection capacity of oncogenic high-risk papillomaviruses (HPV-HR) by multiplex PCR in real time compared to standard swab, taken by a gynecologist or a nurse. With the V-Veil UP<sup>2</sup> veil it is possible to detect oncogenic HPV-HR more than 12 times in 1 woman out of 8.
- ✓ **V-Veil UP<sup>2</sup>**™ is simple, safe, and respects the privacy of all women.

**V-Veil-UP Production** manufactures the best patented self-sampling veil on the market today and the only one to have so many scientific studies which compare the classic method (involving a doctor) with the self-sampling method.

- ✓ **V-Veil UP<sup>2</sup>**™ Self-Sampling Veils are CLASS I Medical Devices
- ✓ Registered in the European database of medical devices (EUDAMED) and the National Agency of Medicines and Medical Devices of Romania (ANMDMR) as a producer of Medical Devices, Systems and Kits.
- ✓ The Self-Sampling Veil is produced in our EN ISO 13485:2016 certified facility and manufacture is done in ISO 8 clean room.
- ✓ Our product is being used WORLDWIDE by a PCR machine manufacturer.
- ✓ **V-Veil UP<sup>2</sup>**™ devices have multiple patent applications and patents pending on several continents.

### SUITABLE FOR:

- 🏠 Doctors    ⚕️ Clinics    🏢 ONGs
- 🚩 National Campaigns
- 🧪 Laboratories    🏥 Hospitals



v-veil-up-production.com

### V-VEIL-UP PRODUCTION

Ciocănești 280B, Călinești  
117197, Argeș, ROMÂNIA

☎ +40 738 886 499

☎ +40 730 108 555

✉ contact@vvup.org



hpv-veil.com

**The Self-Sampling Vaginal Veil V-Veil UP<sup>2</sup>**™ is a class I, non-sterile, UVC disinfected medical device for human use used to collect secretions, various fluids and substances, which are in the vaginal cavity and cervix for diagnosis by laboratory tests.

## COMMERCIAL ADVANTAGES:

- ✓ Suitable for all national HPV detection campaigns.
- ✓ Very practical for NGOs that organize HPV screening campaigns, especially in Muslim or hard to reach communities.
- ✓ Suitable for laboratories that want to have the best product and save money.
- ✓ The only self-sampling veil with scientific studies that show good results, with or without medical assistance.
- ✓ The use of this veil makes it possible to almost double the detection capacity of oncogenic high-risk papillomaviruses (HPV-HR) by multiplex PCR in real time compared to a standard swab method.
- ✓ **V-Veil UP2™** is simple, safe, and respects the privacy of all women.
- ✓ The patented technology **V-Veil UP2™** with its primary tube with **UP2 Retrofitter™**, is suitable for all clinics, hospital and laboratories, as it offers:
  - **Efficiency**
  - **Neutrality with respect to molecular biology reagents**
  - **Compatibility with all HPV PCR Kits on the market**

**V-Veil UP2 14mm**  
Self-Sampling Veil



14mm

Our products are available in:

- individual packaging;
- collective packaging of 250 pcs.;
- bulk at 1000 pcs.

**V-Veil UP2 28mm**  
Self-Sampling Veil



28mm

The 28mm veil can be used to keep medicine/drugs in the vaginal cavity, for a higher efficiency and minimise leaking accidents.

Buy KIT 4.1



Discover all our  
KITS



Discover all our  
scientific studies



**V-Veil UP2 KIT 1.1**  
14mm Self-Sampling Veil



- Content**
- Vaginal Veil Collector V-Veil UP2
  - Instructions for use
  - One container for the storage and transport of samples

**V-Veil UP2 KIT 4.1**  
14mm Self-Sampling Veil



### Content

- Vaginal Veil Collector V-Veil UP2
- Instructions for use
- 1 container for the storage and transport of samples
- 1 safety bag UN3373
- 1 courier bag.


## UP2 RETROFITTER™

### GROUNDBREAKING WORLDWIDE MEDICAL DEVICE INNOVATION

Reduced risk of pre-analytical error and ease of use for automation

The self-sampling medical device **V-Veil UP2™** can be used with its patented, next generation primary tube, **UP2 Retrofitter™**, which allows the laboratory to use this tube with the collected sample, directly in the automation, without the need for an intermediate tube and barcode. **UP2 Retrofitter™** creates free empty space, in which the pipette of the robot can be easily introduced, thanks to the same length and internal diameter of standard hemolysis tube.

Size compatible  
with frozen  
storage boxes



Patent pending

### LABORATORY ADVANTAGES:

- ✓ Safer
- ✓ Faster
- ✓ No contamination
- ✓ Economic
- ✓ Ready for extraction